


VECTOR® TRUEVIEW® AUTOFLUORESCENCE QUENCHING KIT

SKU: SP-8400-15



Description

The TrueVIEW Autofluorescence Quenching Kit provides a novel way to diminish unwanted autofluorescence from non-lipofuscin sources and dramatically improve signal-to-noise ratio. TrueVIEW removes unwanted autofluorescence in tissue sections due to aldehyde fixation, the presence of red blood cells, and structural elements such as collagen and elastin. The quenching action of the kit reagents provides a clear, unambiguous, “true view” localization of the target antigen.

This kit offers you:

- Effectivity with problematic tissues such as kidney and spleen
- An easy-to-use, one-step method
- Compatibility with a wide selection of fluorophores
- Compatibility with standard epifluorescence and confocal laser microscopes
- Availability [with DAPI counterstain](#) and without.

One kit is sufficient to treat approximately 100 to 150 tissue sections.

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.

Specifications


Unit Size	15 ml
Applications	Immunofluorescence, In situ hybridization
Blocking Action	Autofluorescence
Format	Concentrate

Kit Contents:

What's in the box?

- TrueVIEW Reagent A, 5 ml
 - TrueVIEW Reagent B, 5 ml
 - TrueVIEW Reagent C, 5 ml
 - VECTASHIELD Vibrance® Antifade Mounting Medium, 2 ml
-

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.